



WWW.BSMDENTAL.COM

Chengdu Besmile Medical Technology Co., Ltd.

www.bsmdental.com T:+86-28-85317108 E:info@cdbesmile.com Besmile Dental America Inc.

20311 Valley Blvd, Suite# I, Walnut, CA 91789 T: (626)921-5798 E:admin@bsmdentalus.com











Aconia

**INSTRUCTION MANUAL** 

# **Aconia Processing Guide**



**Product Selection** 



CAD



CAM



Coloring



Sintering



Finishing

Aconia Blanks......01
Indication Guide......03
Coloring Liquid......05
General Information...23

Designing......07

 Fast Coloring Technique...........15
Aesthetics Coloring Technique...17


Sintering Parameters...19

Grinding......21
Healing Treatment...22
Sand Blasting......22

# **Product Selection**







**BSM**/Aconia<sup>®</sup>

# Indication Guide

		Veneer	Inlay & Onlay	Reduced crown	Full contour crown	C	Coping	Full contour anterior bridge (3unit)	Full contour posterior bridge (3unit)	Full contour bridge (≤7unit)	Full contour bridge (≤14unit)	Abutment
	TT (3D Multilayer & White)							-86				
НОТ	SHT-ML (3D Multilayer)							•66	<b>W</b> 01	$\wedge$	$\wedge$	
	SHT (White)							•66	<b>W</b> 01	$\wedge$	$\wedge$	
UPGRADED	ST (3D Multilayer & Preshade & White)							•66	<b>W</b> 01	$\wedge$	$\wedge$	
	HT+ (White)							•66	<b>60</b> 0	$\wedge$	$\wedge$	

Ideal results depend on the individual circumstance.

**BSM**/Aconia<sup>°</sup> 04

### **Standard & Master**

### **Artist**



Applied to Aconia® white , TT, SHT ,ST , HT+,specifically.

50ml

50ml



Applied to Aconia® white , SHT ,ST , HT+,specifically.

#### Standard

- Precisely match to VITA\* 16 shade guide
- Dipping & Brushing
- Time saving & Pontic lighting
- · Ideal results applied on Aconia white blanks by Aconia Coloring Technology

The terms marked with \* are registered trademarks and/or brand names of the respective companies.



#### Master

- Precisely matches VITA\* 26 color system
- Dipping & Brushing
- Time saving & Pontic lighting
- · Ideal results applied on Aconia white blanks by Aconia Coloring Technology

The terms marked with \* are registered trademarks and/or brand names of the respective companies.





20ml

#### **Artist**

- Art effect creation
- No need to mix just start
- Individual customization
- · Ideal results applied on Aconia blanks



### Special colors



# Designing

#### Please follow the instructions below:


- 1. To avoid stress concentration, pay attention to path of insertion, undercut, the distance between occlusal and gingival, clearly visible gingival margin.
- 2. Ensure minimum 0.8 mm of labial wall thickness. minimum 1.0 mm of lingual wall thickness.
- 3. Length of continuous pontics under **30** mm Length of continuous cantilevers under 15 mm
- 4. Anterior connector cross section area **9** mm<sup>2</sup> Posterior connector cross section area 12 mm<sup>2</sup>





# **TIPS ON TEETH PREPERATION**

A. Anterior



- Incisal reduction ≥ 2.0 mm
- ② Lingual wall thickness ≥ 1.0 mm
- ③ Labial wall thickness ≥ 0.8 mm
- 4) Outline:

Nature-like Lingual outline Convergence angle 3° to 5°

**B.** Posterior



- Occlusal thickness ≥ 1.3 mm
- ② Occlusal opening angle 120°to140°with rounded edges
- ③ Lingual wall thickness ≥ 1.0 mm
- ④ Labial wall thickness ≥ 0.8 mm
- (5) Convergence angle 6° to 8°









② Cervical margin ≥ 1 mm 90°shoulder with rounded edge

**BSM**/Aconia<sup>®</sup>

### Please follow the instruction of the CAM software

# Milling









#### A. Connector

- ① Set the Connectors on Labial and Lingual sides
- ② Set the connectors on the height of contour between middle 1/3 with cervical 1/3 of the restoration
- ③ No Connectors on the pontic
- (4) No Connectors on the adjacent area

### **B. Sintering frames**

Sintering frames is recommended for the bridge & crowns With more than 7 units



### Loading on the holder

- A. Make sure the holder is clean without any dust before loading
- B. Put the zirconia blank into the holder with the surface horizontal
- C. Gradually tighten the screw in the order 1-2-3-4
- D. Please adjust the position of markers in blank prior to loading, and try to keep the same position each time.

Different Brand/type of milling machines have different procedures in loading, please refer to manufacturer's instructions strictly.

### **De-Spruing**


- A. Professional tool under 10,000-12,000 r/min
- B. Order
- ① Position crowns and bridges on their occlusal surfaces
- ② Start from the connectors on the margin
- 3 End with the connectors on mesial and distal surface

- C. Reduction less than 0.5mm every time
- D. For bridge & Crowns, cut down the connectors outside only
- E. Clean the shaped restoration by brush or air gun.

### Inspection


Cracks or fractures

Faulty restorations must not be processed further.




# **Aconia 3D Multilayer Technology**


### Simply mill, sinter and glaze!




- Gradient chroma:
   Increasing chroma from the top to bottom
- Gradient translucency:
   Increasing translucency from bottom to top
- Gradient flexural strength:
   Increasing flexural strength from top to bottom



Gradient chroma
From the top to bottom



Gradient translucency From bottom to top




Gradient flexural strength


Decrease the tooth wear coefficient

BSM∕Aconia°


# **Aconia 3D Multilayer Technology**




Layer heights:								
14mm	16mm	18mm	20mm	22mm				
2.8 mm	3.2mm	3.6 mm	4 mm	4.4 mm				
2.8 mm	3.2mm	3.6 mm	4 mm	4.4 mm				
7 mm	8mm	9 mm	10 mm	11 mm				
1.4 mm	1.6mm	1.8 mm	2.0 mm	2.2 mm				



Increased body layer - 50% Maximized flexibility of the height usage



Upper the position Higher translucency for the incisal/cusps Occlusal surface parallel to the blank surface



Upper the position Higher translucency for the incisal/cusps Adjust the bridge in the blank to get the ideal result.

# **Fast Coloring Technique**



Methods to avoid pontics discoloration:

Method 1. Decrease the immersing time to 5 sec.

Method 2. Brush diluent before immersing.



Immersing time can be modified according to the wall thickness of the restorations.

# **Aesthetics Coloring Technique**

























Shake well

**Aconia Artist** Aconia Standard/Master prepared ①Brush TO1

② Brush TO2/TO3

(3) Immersing

4 Dabbing

⑤ Drying

6 Brush O1&O2 & Special colors

7 Drying

On Incisal/Cusp X 2-3


10-15 sec. for single crown 5-8 sec. for bridge

10min at 80 - 100°C

30min for single crown 60min for bridge at 80 - 100°C







The coloring skills will be enhanced with the experience accumulation of dental technicians.

**BSM**/Aconia<sup>®</sup>

17

### 

1520°C

370

(4.4°C/min)

90min

1520°C

460

(-3°C/min)

(MIN.)

800°C

### Aconia Sintering Parameters (1-7 units)

Applied to ST/SHT/TT(Multilayer), ST(Preshade), HT+/ST/SHT(White)

Step	Initial Temp.	Rate	Final Temp.	Time		
1	50°C	2.8°C/min	300°°C	90min		
2	300℃	5°C/min	1520℃	240min		
3	1520℃	Holding	1520℃	90min		
4	1520℃	- 4.5°C/min	800°C	160min		
5	800	Natural cooling				


<sup>\*</sup>Applied to single unit and bridges under 7 units(1-7units)

### Aconia Bridge Sintering Parameters(>7 units)

Applied to ST/SHT(Multilayer), ST(Preshade), HT+/ST/SHT(White)

Step	Initial Temp.	Rate	Final Temp.	Time		
1	50°C	2.8°C/min	300°C	90min		
2	300°C	4.5°C/min	1520°C	280min		
3	1520°C	Holding	1520℃	90min		
4	1520°C	- 3°C/min	800°C	240min		
5	800	Natural cooling				

<sup>\*</sup> Applied to bridges above 7 units(>7units)



1530°C


55

30min

1530°C

85





\*It's recommended to use BSM-S30 sintering furnace.Input "-121" to enter into the natural cooling process.

# Aconia® Fast Sintering Parameters (1-3 units)

Applied to **ST/SHT/TT**(Multilayer),**ST**(Preshade), **HT+**(White)

Step	Initial Temp.	Rate	Final Temp.	Time		
1	40°C	34°C/min	900°C	25min		
2	900°C	21°C/min	1530°C	30min		
3	1530℃	Holding	1530℃	30min		
4	1530℃	- 27°C/min	1000°C	20min		
5	1000°C	Natural cooling				

<sup>\*</sup>Applied to single unit and bridges under 3 units(1-3 units)

\*It's recommended to use BSM-FC30 sintering furnace.Input"-20"to enter into natural cooling process.

B. It has been proved that above curves lead to the best outcome of Aconia®, and sintering programs also can be adjusted according to the specific situations, and different furnaces.

(MIN.)

105

C. Please clean the furnace regularly.

25

(21°C/min)

(34°C/min)







<sup>\*</sup>It's recommended to use BSM-S30 sintering furnace.Input "-121" to enter into the natural cooling process.

<sup>\*</sup>It's recommended to use BSM-S30 sintering furnace.Input "-121" to enter into the natural cooling

<sup>\*</sup>For the full-arch restorations/ thick-crown restorations etc., it's suggested to extend the holding

<sup>\*</sup>Long span bridges should be sintered with sintering frame

A. Sintering can take place in all common dental sintering furnaces which are approved for the sintering of restorations made of zirconium dioxide.

**Finishing** 

# **Finishing**



### Grinding is not recommended after sintering, if it is inevitable, please follow the tips as below:

(1) Professional dental zirconia grinding instrument




















- (3) Grinding with a light pressure by a single direction
  - (4) Grinding speed 20,000-25,000 r/m
  - (5) No grinding inside the crown

### **Healing treatment**

Recommend to increase the strength and blinding force to zirconium porcelain



#### Sand blasting

Sand blasting is recommended for cleaning, recovering material strength as well as a better connection with veneering material.

- 1. Before sand blasting, cleaning is necessary with ultrasound in a water bath or a steam jet.
- 2. High purity alumina sand —50 µ m(270 mesh)
- 3. Under the pressure of 2-2.5Bar(0.2-0.25mpa)

For veneering the restoration, staining and glazing, as well as cut-back and layering techniques,

or a combination of both, are suitable for Aconia

Please follow the instruction for use of the manufacturer.

**BSM**/Aconia®

# **Technical Data**

	TT-ML	SHT-ML	ST-ML	TΤ	ST (Preshade)	HT+
Translucency	43%-49%	43%-46%	42%-45%	49%	43%	42%
Flexural strength/Mpa	600-900	900-1100	1000-1200	800±100	1250±100	1400±100
Density/(g/cm³)	>3	>3	>3	>3	>3	>3
Sintered Density/ (g/cm³)	>6.02	>6.02	>6.02	>6.02	>6.02	>6.02
Fracture toughness/ (Mpa.m <sup>1/2</sup> )	>3	>5	>5.5	>3	>5.5	>9
Vickers- hardness HV10	1300±5%	1300±5%	1300±50	1300±5%	1300±5%	1300±5%
СТЕ	(10.5±0.5)* 10 <sup>-6</sup> K <sup>-1</sup>					
Radioactivity/Bq·g <sup>-1</sup>	<0.1	<0.1	<0.1	<0.1	<0.1	<0.1
Chemical solubility/ (µg/cm²)	<50	<50	<50	<50	<50	<50

# Storage & Handling

Please check the package and product carefully before using. Operated only by professionals.

Transportation avoid damage from collision, extrusion and vibration.

Storage: The product must be stored in its original packaging at room temperature.

Expiry date: Long-term.

### Warning

Operator must wear dust mask.

Do not get dust into eyes or directly contact with skin.

### Symbol:

M	Date of Production	X	Expiry date	LOT	Batch number
$\triangle$	Caution		Consult instruction for use	Ţ	Fragile
T T	Keep dry				

Registrant, Production enterprise name: Chengdu Besmile Biotechnology Co., Ltd.

After sales service by: Chengdu Besmile Biotechnology Co., Ltd.