


■ Global Reach, Local Precision


Business Layout: With subsidiaries in the USA and Europe, BSM serves clinicians and labs in over 100 countries and regions.



Certifications: Our commitment to quality is validated by NMPA (China), FDA (USA), CE (Europe), ISO 13485, and numerous other international certifications.



R&D Focus: Significant investment (>20% annual revenue) drives innovation in materials science and precision engineering, supported by a dedicated R&D team comprising ≈20% of our 400+ staff.



About Besmile

Besmile is a global digital solutions provider for restorative and implant dentistry.

We offer integrated systems—premium CAD/CAM materials, advanced equipment, and precision implant systems—designed to streamline workflows and ensure reliable results.

All core products are developed and manufactured in-house, ensuring consistent quality and continuous innovation.

Trusted by over 1,000 partners in more than 100 countries, we empower dental professionals to create confident, lasting smiles.

Technology creates the best smile.



Chengdu Besmile Biotechnology Co.,Ltd.

www.bsmdental.com

T: +86-28-85317108

E: info@cdbesmile.com

Besmile Dental America Inc.

20311 Valley Blvd, Suite#1, Walnut, CA91789

T: (626)921-5798

E: sales@bsmdentalus.com



Contact Us



BSM Dental Implant Components

Your Source for Innovation & Quality



Partnering for Clinical Excellence

Welcome to the BSM Dental Implant Components . BSM Dental is dedicated to providing superior implant restorative solutions through advanced R&D and precision manufacturing.

For our distributor partners, we offer a market-competitive product line to drive business growth. For dental laboratories, we provide high-precision, easily integrated components that optimize and streamline restorative workflows, ensuring outstanding clinical results.

Our commitment is to support dental professionals in achieving predictable, aesthetic, and long-lasting patient outcomes. Explore our comprehensive range designed for your clinical success.



Elevating Efficiency & Expanding Opportunities: The BSM Advantage

Wide Range of Compatibility:

Seamless integration with >40 major implant brands and >100 interface types.



Triangle a



Triangle b



Crossfit



SynOcta



Hexagonal



Cone











BSM Dental Implant Components

Precision Compatibility for Predictable Restorations



Full range of Implant Components

Full range of standard components, custom abutment services, and digital workflow solutions (scanbodies, libraries) are included.




Premium Materials:

Utilizing Titanium Grade 5 ELI (Ti-6Al-4V ELI) sourced from globally leading suppliers (e.g., Carpenter Technology), ensuring strength (>1000Mpa) and biocompatibility.



Exceptional Precision:

- Fit clearance < 10µm, reducing micro-leakage risks.
- Achieved through multi-axis CNC machining with micron-level accuracy.



Rigorous Quality Control:

- 8 distinct inspection processes from raw material to finished product.
- 100% interface inspection guarantees perfect fit and profile matching.



High Quality:

- >5 Million abutments successfully placed worldwide.
- Clinical data shows incidence of technical complications < 0.02%, significantly below industry averages.